
















Lab Driver




- Convenient for use with untouchable driver part
- Easy to exchange tips
- Heat treated
- 'Brinell hardness' is higher than 55



** Two different LAB Driver lengths of 30mm*

SIZE	SHAPE	COLOR	RECOMMENDED TORQUE
 Ball Torx		Orange	20 Ncm
 Ball Hex		White	20 Ncm
 Hex 1.0		Blue	15 Ncm
 Hex 1.2		Purple	20 Ncm
 Hex 1.25		Black	30 Ncm
 Hex 1.27		Red	30 Ncm
 Hex 1.3		Green	30 Ncm
 Torx 1.5		Yellow	35 Ncm
 Torx 1.7		Brown	35 Ncm
 Cross		Gray	30Ncm

Clinic Ball Driver



- Ball Driver Compatible with Straumann® Ratchet Wrench
- Available in surgical stainless steel
- Can be cleaned easily and thoroughly
- Heat treated
- For angled screw channel
- 'Brinell hardness' is higher than 55

SIZE	LENGTH	ORDER CODE	IMPLANT SYSTEM	RECOMMENDED TORQUE
 Ball Torx	22mm	BT1722	1. Most of angled screws of ARUM Ti-base (screw head size is more than 2.1mm) 2. All implant system for disk angulated screw channels	20 Ncm
	15mm	BT1715		
	8mm	BT1708		
 Ball Hex	22mm	BH1222	Some of angled screws of ARUM Ti-base (screw head size is less than 2.1mm) *Please contact us for a separate list.	20 Ncm
	15mm	BH1215		
	8mm	BH1208		

