

ARUM Multi-Unit

ARUM Product List Multi-Unit



No	Arum® Product	Compatible Product Group	Pg	Screw	Analog	Ti-Base	Model/ Intraoral Scanbody	Pre-milled Blank 10mm
1	AD-MT	ADIN®	140	○	-	-	○	-
2	AN-MT	ANTHOGR® axiom®	140	○	-	-	○	-
3	AT-MT	ASTRA TECH™	141	○	-	-	○	-
4	AT-MT	ASTRA TECH™ OsseoSpeed™ EV	141	○	-	-	○	-
5	BEG-MT	BEGO®	142	○	-	-	○	-
6	BM3I-MT	BIOMET 3i®	142	○	-	-	○	-
7	BH-MT	BIOHORIZONS®	143	○	○	-	○	-
8	BRM-MT	BRENT Medical SKY®	143	○	-	-	○	-
9	BT-MT	BTI	144	○	○	-	○	-
10	BT-MT	BTI Angled converter	144	○	-	-	○	-
11	CON-MT	CONELOG®	145	○	○	-	○	-
12	CRT-MT	CORTEX™	145	○	-	-	○	-
13	DTS-MT	DENTIS®	146	○	-	-	○	○
14	DTM-MT	DENTIUM®	147	○	○	○	○	○
15	DS-ANK	DENTSPLY® Ankylos®	148	○	○	○	○	○
16	DS-MT	DENTSPLY® Xive®	149	△	△	-	△	-
17	DO-MT	DIO® UF	149	○	○	-	○	-
18	EU-MT	EUROTEKNIKA™	150	△	-	-	○	-
19	GLD-MT	Global D®	150	○	-	-	○	-
20	ID-MT	IMPLANT DIRECT®	151	○	-	-	○	-
21	JD-MT	JD EVOLUTION®	151	○	-	-	○	-
22	MM-MT	Medentis Medical® ICX	152	○	-	-	○	-
23	MS-MT	MIS®	152	○	-	-	○	-
24	MG-MT	MEGAGEN®	153	△	-	-	○	-
25	MG-MT	MEGAGEN® Octa	153	△	-	-	○	-

ARUM Product List Multi-Unit



No	Arum® Product	Compatible Product Group	Pg	Screw	Analog	Ti-Base	Model/ Intraoral Scanbody	Pre-milled Blank 10mm
26	MG-MT	MEGAGEN® N-Type	154	○	-	-	○	-
27	NED-MT	NEODENT®	154	○	-	-	○	-
28	NB-MT	NOBEL BIOCARE®	155	○	○	○	○	-
29	OST-MT	OSSTEM®	157	○	-	-	○	-
30	SC-MT	SIC	158	○	-	-	○	-
31	SI-MT	SOUTHERN IMPLANTS	158	○	-	-	○	-
32	STR-MT	STRAUMANN® Bone Level®	159	○	○	-	○	-
33	STR-MT	STRAUMANN® Bone Level® Angle	159	○	-	-	○	-
34	STR-MT	STRAUMANN® synOcta®	160	○	-	-	○	-
35	STR-SRA	STRAUMANN® Screw Retained Abument®	161	○	○	○	○	○
36	ZMR-TSV	ZIMMER® Tapered Screw-Vent®	160	○	-	-	○	-




△ will be developed within 6 months, please contact us for further details.

compatible with

ADIN®



MULTI-UNIT

	Screw	MS244 [HEX1.27]
	Screw <i>Pack of 10</i>	PMS244
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M218_02
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM218_02

compatible with




ANTHOGRYR® axiom®



AXIOM-MULTI-ANGLED 4.8



AXIOM-CONVERTIBLE 4.0

	Screw	MS222 [HEX 1.2]	MS223 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS222	PMS223
	Recommended torque	15Ncm	15Ncm
	Model/Intraoral Scanbody	M222	M223
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM222	PM223

compatible with

ASTRA TECH™



ASTRATEC-UNI 2.0



ASTRATEC-UNI 3.0



ASTRATEC-UNI 4.5




	Screw	MS209 [HEX 1.27]	MS672 [HEX 1.3]	MS209 [HEX 1.27]
	Screw <i>Pack of 10</i>	PMS209	PMS672	PMS209
	Recommended torque	15Ncm	15Ncm	15Ncm
	Model/Intraoral Scanbody	M209	M242	M210
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM209	PM242	PM210

compatible with

ASTRA TECH™








ASTRATEC ANGLED(TCA)

	Screw	MS209 [HEX 1.27]
	Screw <i>Pack of 10</i>	PMS209
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M211
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM211

compatible with

BEGO®






MULTI-UNIT

	Screw	MS208 [HEX 1.27]
	Screw <i>Pack of 10</i>	PMS208
	Recommended torque	20Ncm
	Model/Intraoral Scanbody	M208
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM208

compatible with

BIOMET 3i®



BIOMET-MULTI

	Screw	MS208 [HEX 1.27]
	Screw <i>Pack of 10</i>	PMS218
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M218
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM218

compatible with

BIOHORIZONS®

MULTI-UNIT 4.8

	Screw	MS245 [HEX1.25]
	Screw <i>Pack of 10</i>	PMS245
	Recommended torque	15Ncm
	Analog	MF200
	Analog <i>Pack of 10</i>	PMF200
	Model/Intraoral Scanbody	M200_02
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM200_02






compatible with

**BREDENT Medical
SKY®**

SKY-UNI.CONE



SKY-FAST&FIXED

	Screw	MS229 [TORX]	MS229 [TORX]
	Screw <i>Pack of 10</i>	PMS229	PMS229
	Recommended torque	18Ncm	18Ncm
	Model/Intraoral Scanbody	M228	M229
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM228	PM229





compatible with
BTI



Multi-Unit 4.0



Multi-Unit 5.5

	Screw	MS241 [HEX 1.2]	MS241 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS241	PMS241
	Recommended torque	20Ncm	20Ncm
	Analog	MF241	MF241
	Analog <i>Pack of 10</i>	PMF241	PMF241
	Model/Intraoral Scanbody	M241	M241
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM241	PM241

compatible with
BTI
Angled Converter



5.0

	Screw	MS241 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS241
	Recommended torque	20Ncm
	Model/Intraoral Scanbody	M221
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM221

compatible with
CAMLOG®
CONELOG®



MULTI-UNIT 3.3



MULTI-UNIT 3.8



MULTI-UNIT 4.3






MULTI-UNIT 5.0

	Screw	MS251	MS251	MS251	MS252
	Screw <i>Pack of 10</i>	PMS251	PMS251	PMS251	PMS252
	Recommended torque	15Ncm	15Ncm	15Ncm	15Ncm
	Analog	MF242	MF242	MF242	MF243
	Analog <i>Pack of 10</i>	PMF242	PMF242	PMF242	PMF243
	Model/Intraoral Scanbody	M251	M251	M251	M252
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM251	PM251	PM251	PM252

compatible with
CORTEX™







CORTEX-MULTI

	Screw	MS232 [HEX 1.27]
	Screw <i>Pack of 10</i>	PMS232
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M232
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM232

compatible with

DENTIS®

Multi unit Angled

	Screw	MS234 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS234
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M234
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM234
	Pre-milled Blank 10mm	MA234









compatible with

DENTIUM®

CONVERTIBLE 4.5



CONVERTIBLE 5.5

	Screw	MS212 [HEX 1.27]	MS212 [HEX 1.27]
	Screw <i>Pack of 10</i>	PMS212	PMS212
	Recommended torque	15Ncm	15Ncm
	Analog	MF212	MF239
	Analog <i>Pack of 10</i>	PMF212	PMF239
	Ti-base	MTi_212	-
	Angled screw driver type	Torx	-
	Type of sleeve for casting and pressing	F	-
	Model/Intraoral Scanbody	M212	M239
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM212	PM239
	Pre-milled Blank 10mm	MA212	-









compatible with

DENTSPLY®
Ankylos®

BALANCE BASE NARROW



BALANCE BASE

	Screw	MS235 [HEX 1.0]	MS236 [HEX 1.0]
	Screw <i>Pack of 10</i>	PMS235	PMS236
	Recommended torque	15Ncm	15Ncm
	Analog	MF235	MF236
	Analog <i>Pack of 10</i>	PMF235	PMF236
	Ti-base	MTi_235	MTi_236
	Angled screw driver type	Torx	Torx
	Type of sleeve for casting and pressing	B	C
	Model/Intraoral Scanbody	M235	M236
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM235	PM236
	Pre-milled Blank 10mm	MA235	MA236

compatible with

DENTSPLY®
Xive®

3.4







3.8



4.5







5.5

	Screw	MS238 [HEX1.22]	MS239 [HEX1.22]	MS240 [HEX1.22]	MS248 [HEX1.22]
	Screw <i>Pack of 10</i>	PMS235	PMS239	PMS240	PMS248
	Recommended torque	24Ncm	24Ncm	24Ncm	24Ncm
	Analog	MF244	MF245	MF246	MF247
	Analog <i>Pack of 10</i>	PMF244	PMF245	PMF246	PMF247
	Model/Intraoral Scanbody	M243	M244	M245	M246
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM243	PM244	PM245	PM246

compatible with

DIO®
UF






MULTI-UNIT 4.8

	Screw	MS227 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS227
	Recommended torque	20Ncm
	Analog	MF227
	Analog <i>Pack of 10</i>	PMF227
	Model/Intraoral Scanbody	M227
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM227

compatible with

EUROTEKNIKA™






MULTI-UNIT 4.8

	Screw	M243
	Screw <i>Pack of 10</i>	PM243
	Recommended torque	20Ncm
	Model/Intraoral Scanbody	M218_03
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM218_03

compatible with

IMPLANT DIRECT®

MULTI-UNIT 5.0

	Screw	MS214 [HEX 1.27]
	Screw <i>Pack of 10</i>	PMS214
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M214
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM214






compatible with

GLOBAL D®

MULTI-UNIT SYSTEM STRAIGHT








MULTI-UNIT SYSTEM-ANGLED

	Screw	MS215 [HEX 1.24]	MS216 [HEX 1.24]
	Screw <i>Pack of 10</i>	PMS215	PMS216
	Recommended torque	15Ncm	15Ncm
	Model/Intraoral Scanbody	B215	B216
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM215	PM216

compatible with

**J DENTAL CARE®
JD EVOLUTION®**






MULTI-UNIT

	Screw	MS247 [HEX1.2]
	Screw <i>Pack of 10</i>	PMS247
	Recommended torque	20Ncm
	Model/Intraoral Scanbody	M200_04
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM200_04

compatible with
Medentis Medical®
ICX








ICX-MULTI

	Screw	MS233 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS233
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M200_01
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM200_01

compatible with
MIS®








MULTI-UNIT SYSTEM STANDARD

	Screw	MS217 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS217
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M217
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM217

compatible with
MEGAGEN®








MULTI

	Screw	MS250 [HEX1.2]
	Screw <i>Pack of 10</i>	PMS250
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M221_01
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM221_01

compatible with
MEGAGEN®
Octa






MULTI 4.8

	Screw	MS249 [HEX1.2]
	Screw <i>Pack of 10</i>	PMS249
	Recommended torque	30Ncm
	Model/Intraoral Scanbody	M202_01
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM202_01

compatible with
MEGAGEN®
N-Type



MULTI 4.8

	Screw	MS246 [HEX1.2]
	Screw <i>Pack of 10</i>	PMS246
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M200_03
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM200_03

compatible with
NEODENT®
CM



MICRO Conical



MINI Conical

	Screw	MS243 [HEX1.2]	MS243 [HEX1.2]
	Screw <i>Pack of 10</i>	PMS243	PMS243
	Recommended torque	10Ncm	10Ncm
	Model/Intraoral Scanbody	M217-NED	M218_00-NED
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM217-NED	PM218_00-NED




compatible with
NEODENT®
GM



MICRO Conical



MINI Conical

	Screw	MS243 [HEX1.2]	MS243 [HEX1.2]
	Screw <i>Pack of 10</i>	PMS243	PMS243
	Recommended torque	10Ncm	10Ncm
	Model/Intraoral Scanbody	M217-GM	M218_00-GM
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM217-GM	PM218_00-GM








compatible with
NOBEL BIO CARE®



MULTI-UNIT 4.8










MULTI-UNIT 6.0

	Screw	MS200 [TORX]	MS201 [TORX]
	Screw <i>Pack of 10</i>	PMS200	PMS201
	Recommended torque	15Ncm	15Ncm
	Analog	MF200	MF201
	Analog <i>Pack of 10</i>	PMF200	PMF201
	Ti-base	MTI_200	MTI_201
	Angled screw driver type	Torx	Torx
	Type of sleeve for casting and pressing	E	F
	Model/Intraoral Scanbody	M200	M201
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM200	PM201

compatible with

NOBEL BIOCARE®

MULTI-UNIT PLUS 4.8

	Screw	MS200_01 [TORX]
	Screw <i>Pack of 10</i>	PMS200_01
	Recommended torque	15Ncm
	Analog	MF680
	Analog <i>Pack of 10</i>	PMF680
	Ti-base	MTI_680
	Angled screw driver type	Torx
	Type of sleeve for casting and pressing	E
	Model/Intraoral Scanbody	M218_01
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM218_01




compatible with

OSSTEM®

TS MULTI ANGLED Mini



TS MULTI ANGLED Regular

	Screw	MS213 [HEX 1.2]	MS213 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS213	PMS213
	Recommended torque	15Ncm	15Ncm
	Model/Intraoral Scanbody	M213	M213
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM213	PM213

compatible with

OSSTEM®



TS CONVERTIBLE 4.0



TS CONVERTIBLE 5.0








TS CONVERTIBLE 6.0

	Screw	MS224 [HEX 1.2]	MS225 [HEX 1.2]	MS226 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS224	PMS225	PMS226
	Recommended torque	15Ncm	15Ncm	15Ncm
	Model/Intraoral Scanbody	M224	M225	M226
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM224	PM225	PM226

compatible with

SIC

SIC-MULTI

	Screw	MS221 [HEX 1.2]
	Screw <i>Pack of 10</i>	PMS221
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M221
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM221






compatible with

SOUTHERN IMPLANTS®

MULTI-UNIT 4.8










MULTI-UNIT 6.0

	Screw	MS242 [HEX1.22]	MS242 [HEX1.22]
	Screw <i>Pack of 10</i>	PMS242	PMS242
	Recommended torque	10-15Ncm	10-15Ncm
	Model/Intraoral Scanbody	M200_05	M201_01
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM200_05	PM201_01

compatible with

STRAUMANN®
Bone Level®NC_CONVERTIBLE
D3.5NC_CONVERTIBLE
D4.5RC_CONVERTIBLE
D4.5RC_CONVERTIBLE
D6.5

	Screw	MS204 [TORX]	MS205 [TORX]	MS205 [TORX]	MS206 [TORX]
	Screw <i>Pack of 10</i>	PMS204	PMS205	PMS205	PMS206
	Recommended torque	15Ncm	15Ncm	15Ncm	15Ncm
	Analog	MF204	MF205_NC	MF205_RC	MF206
	Analog <i>Pack of 10</i>	PMF204	PMF205_NC	PMF205_RC	PMF206
	Model/Intraoral Scanbody	M204	M205	M205	M206
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM204	PM205	PM205	PM206






compatible with

STRAUMANN®
Bone Level® Angle

NC 2.5



RC 2.5

	Screw	MS207 [TORX]	MS207 [TORX]
	Screw <i>Pack of 10</i>	PMS207	PMS207
	Recommended torque	15Ncm	15Ncm
	Model/Intraoral Scanbody	M207	M207
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM207	PM207



compatible with

STRAUMANN®
synOcta®

RN SynOcta






WN SynOcta

	Screw	MS202 [TORX]	MS203 [TORX]
	Screw <i>Pack of 10</i>	PMS202	PMS203
	Recommended torque	15Ncm	15Ncm
	Model/Intraoral Scanbody	M202	M203
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM202	PM203

compatible with

ZIMMER®
Tapered Screw-Vent®

Tapered Abutment

	Screw	MS671 [HEX 1.27]
	Screw <i>Pack of 10</i>	PMS671
	Recommended torque	15Ncm
	Model/Intraoral Scanbody	M240
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM240




compatible with

STRAUMANN®
Screw Retained Abutment®

Ø3.5mm



Ø4.6mm

	Screw	MS237 [TORX]	MS238 [TORX]
	Screw <i>Pack of 10</i>	PMS237	PMS238
	Recommended torque	15Ncm	15Ncm
	Analog	MF237	MF238
	Analog <i>Pack of 10</i>	PMF237	PMF238
	Ti-base	MTi_237	MTi_238
	Angled screw driver type	Torx	Torx
	Type of sleeve for casting and pressing	C	E
	Model/Intraoral Scanbody	M237	M238
	Model/Intraoral Scanbody <i>Pack of 5</i>	PM237	PM238
	Pre-milled Blank 10mm	MA237	MA238